

Alternative Drive Controls: Giving Power to the Masses

Presented by Jacelyn Goh

Clinical Educator, OT and Product Consultant
Linds Rehabilitation Equipment
jacelyn@lindsrehab.com.au



Orlando

Learning Objectives

- Participants will be able to list methods to improve driving skills through standard set up prior to considering alternative driving systems
- Participants will be able to differentiate different driving systems for powered wheelchairs and list the benefits for different client groups
- Participants will be to list 2 considerations for additional functions of alternative driving systems

Alternate Drive Controls – What and Why?

Standard Drive System

- R-Net and Q-Logic
- Proportional control
- With standard joystick/modified joystick
- Mounted on armrest/other suitable locations
- Expandable vs Non-Expandable
- Function controls on the controller/Easy access

Physical requirements:

- Sufficient UL movement and fine motor control
- Ability to sustain force
- Joystick placement




Q-LOGIC 3
Advanced Drive Control System



Alternate Drive Controls – What and Why?

Adjustability within this standard system:

- Programming (Speed, Acceleration, Active Throw, Dead Band, Tremor dampening, Active Orientation)
- Joystick Handles
- Mounting systems
- Access buttons




Alternate Drive Controls – What and Why?

Alternate Drive Systems

- Different sources of control
- Placement of systems
- Would require additional control interface

Why use an alternate drive system?

- To accommodate for client's needs/ability to drive
- Better posture/access to support endurance



<https://www.numotion.com/products-services/innovative-lifestyle-products/ability-drive#overview>



Assessment Process for Alternate Drive Controls:

Physical Abilities

- Available Movements (ULs / Head/ Oral)
- Vision
- Postural Stability and Displacement
- Prognosis and impact on abilities

Wheelchair requirements

- Power Functions
- Drive wheel location
- Communication devices

Access

- Environment use

Cognition / Perception

- Ability to learn
- Good safety awareness
- Problem-solving
- Spatial perception/awareness

Carer Requirement

- Ability to manage complex systems
- Transfers / daily activities in wheelchair

Alternate Control Drive System - Based on Consistent Movement Patterns

Good / Fair Gross Motor Control, Limited Fine Motor

- **Switch (4) / Proximity Switch**
- Single Switch scanning (single body site)
- CoMove It



JSC for foot drive (FCMH)




JSC for arm drive (ACMH)

Alternate Control Drive System - Based on Consistent Movement Patterns

Good/Fair Gross Motor Control, Limited Fine Motor

- Switch (4) / Proximity Switch
- **Single Switch scanning (single body site)**
- CoMove It

<https://hub.permobil.com/blog/alternative-drive-controls-single-switch-scanning-eye-gaze>



<https://www.youtube.com/watch?v=UAudJSUYDDo>



Alternate Control Drive System - Based on Consistent Movement Patterns

Good/Fair Gross Motor Control, Limited Fine Motor

- Switch (4) / Proximity Switch
- Single Switch scanning (single body site)
- **CoMove It**



Alternate Control Drive System - Based on Consistent Movement Patterns

Good Fine motor, Limited Gross motor


- Mini and Micro- Switches
 - Significantly lesser amount of force for displacement (~10-50g)
 - Multiple mounting options



Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- **Chin control (micro switches)**
- Gyro systems (Vigo and Munevo Glasses)
- Head Controls (proportional vs Head array)



Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- Chin control (micro switches)
- **Gyro systems (Vigo and Munevo Glasses)**
- Head Controls (proportional vs Head array)



<https://www.sunrisemedical.com/power-wheelchairs/switch-it-electronics/head-controls/vigo>




Alternate Control Drive System - Based on Consistent Movement Patterns

No UL, Good Head Control

- Chin control (micro switches)
- Gyro systems (Vigo and Munevo Glasses)
- **Head Controls - Switched Head array (Permobil total control, Switch-it Duo Pro, ASL Fusion, Stealth i-Drive)**



SUNRISE MEDICAL Dual Pro™



Stealth Products®
i-Drive® Head Arrays

Alternate Control Drive System - Based on Consistent Movement Patterns

Fair Head Control, Limited Extremity control

- Switched Head Array
- Munevo Glasses

Considerations:

- Head supports required
- Endurance
- Increased adjustments to profiles to reduce postural displacement

Alternate Control Drive System - Based on Consistent Movement Patterns

Little head control, Good oral control

- Sip n Puff, Combination



Alternate Control Drive System - Based on Consistent Movement Patterns

Minimum head control, good eye gaze control

- Eye gaze



Mounting Systems

Setting up the most ergonomic position of drive method to ensure endurance and ease of use

Considerations for carer staff:

- Ease to remove/replace
- Adjustments required
- Durability

Mounting Systems



Tray mounted Systems



Arm rest mounts



Center Drive System/Mounts



<https://activecontrols.com/info/cds>



Chin Control / Head Array mounts



Recommendation

“The factors identified in research related to non-use of provided assistive technology (as listed) should be considered by the therapist during the wheelchair prescription process, as these may influence the outcomes.”

https://www.aci.health.nsw.gov.au/__data/assets/pdf_file/0003/167286/Guidelines-on-Wheelchair-Prescription.pdf Page 25



You made it!

Questions?

Please come and try out the systems

Jace Goh
Occupational Therapist/ Product Consultant
jacelyn@lindsrehab.com.au



Cubro[®]