

# Rapid AMR diagnostics via recombinase polymerase amplification – our experiences to date

Dr Nicole Ertl

The University of Queensland – Centre for Clinical Research

# Background

- Early detection of antimicrobial resistance (AMR):
  - Targeted treatment
- Potential of isothermal amplification tests closer to/at point-of-care
  - Fast turn-around time, sensitive and specific, versatile, easy to carry out and read, no specialised equipment
- Our experiences developing AMR tests using recombinase polymerase amplification (RPA)

# Methods


## Targets

- *Neisseria gonorrhoeae*<sup>1</sup>
- Ciprofloxacin susceptibility<sup>1</sup>
- ESBLs<sup>2</sup>
- AmpC<sup>2</sup>
- Carbapenemases


## Extract DNA

- E.g.
- Commercial kit
  - Boiled preparation
  - Rapid extraction reagent

## Amplify



## Detect



## Time:

**10-30 min**

**10-15 min**

**5 min**

# Results

Test parameters	<i>N. gonorrhoeae</i> <sup>1</sup>	Ciprofloxacin <sup>1</sup>	ESBLs <sup>2</sup>	AmpC <sup>2</sup>	Carbapenemases
Diagnostic target(s)	<i>porA</i>	<i>gyrA</i> wildtype	<i>bla</i> <sub>CTX-M group 1</sub> <i>bla</i> <sub>CTX-M group 9</sub>	<i>bla</i> <sub>CMY-2</sub>	<i>bla</i> <sub>KPC</sub> <i>bla</i> <sub>IMP</sub> <i>bla</i> <sub>OXA-48</sub>
Amplification temp (°C)	39	39	39	39	39
Detection	lateral flow	lateral flow	lateral flow	lateral flow	lateral flow
Analytical sensitivity (copies/reaction)	11.4	247.4	≤ 27.3	3.4	≤ 17.7
Clinical sensitivity (%) <sup>a</sup>	87.5	83.3	100*	100	-
Clinical specificity (%) <sup>a</sup>	-	100	≥ 94.1	100	≥ 97.9
Time to result (min)	15	15	15	15	15-20

<sup>a</sup>compared to PCR and/or Sanger sequencing

\*only *bla*<sub>CTX-M group 1</sub> clinical sensitivity could be determined



- Similar sensitivity and specificity as PCR
- Rapid (sample to result:  $\leq 30$  min)
- Easy to adapt
- Easy readout
- No specialised equipment / hand incubation<sup>1</sup>



- SNP detection difficult (tolerant to mismatches)
- Increase user friendliness by decreasing number of steps
- Manufacturing issues

# Conclusion

- Potential as closer to/at point of care test:
  - Rapid sample to result when coupled with rapid extraction method
  - No specialised equipment or extensive training needed
- Potential to increase testing accessibility in low resource environments and low income countries
- Continue collaborating with industry with aim to take our research into practice

# Disclosure of interest

Work was conducted in conjunction with BioCifer Pty Ltd and supported by a Children's Hospital Foundation grant and funded by DMTC. JM is a shareholder and director of BioCifer.



## Contact

**Dr Nicole Ertl**

[nicole.ertl@uq.edu.au](mailto:nicole.ertl@uq.edu.au)

## Acknowledgements

Dr Jacob Tickner

Dr Adam Irwin

Mr Abdulrahman Ayfan

Dr Brian Forde

Ms Michelle Bauer

Ms Claire Wang

Dr Patrick Harris

Dr Claire Heney

Dr Hosam Zowawi

Prof David Paterson

Prof Monica Lahra

A/Prof Joanne Macdonald

A/Prof David Whiley

Queensland Paediatric Infectious Diseases