

# The Deviron® Detergent Portfolio

Our Emprove® detergents for a sustainable Triton™ X-100 replacement in biologics manufacturing

Biologics manufacturing processes include critical steps such as viral inactivation or cell lysis that are typically performed with detergents.

Triton™ X-100 detergent (4-tert-octylphenol polyethoxylate) is the primary component templated in biomanufacturing processes, but as of January 2021, its unauthorized use has been prohibited in the European Union by the European Commission, due to its listing to REACH<sup>1</sup> (Registration, Evaluation, Authorization and Restriction of Chemicals).

In response to the global implications of restrictions on Triton™ X-100, we have developed the Deviron® C16 Emprove® Evolve detergent and the Deviron® 13-S9 Emprove® Expert detergent.

**A comparison table is available to aid product evaluation and selection.**

	Deviron® C16 Emprove® Evolve Detergent	Deviron® 13-S9 Emprove® Expert Detergent	Triton™ X-100 Emprove® Expert Detergent
Product code	108693	108694	108643
Surfactant type	Zwitterionic (pI 8.9)	Non-ionic	Non-ionic
Chemical name	N,N-Dimethyltetradecylamin-N-oxide	Alcohols, C11-15-secondary, ethoxylated	4-(1,1,3,3-tetramethylbutyl)phenol ethoxylated
CAS number	3332-27-2	68131-40-8	9036-19-5
Critical micelle concentration (CMC)	0.002-0.003% wt (24 °C)	0.005% wt (24 °C)	0.013% wt (24 °C)
Biodegradability (Readily biodegradable: >60% within 28 days) OECD 301	Readily biodegradable: 88% degradation (OECD 301B)	Readily biodegradable: 74% degradation (OECD 301B)	Substance of Very High Concern as defined by REACH <sup>1</sup>
Virus inactivation efficiency LRV ≥ 4 (@ ≥0.5% wt, 15-25 °C)	Better or comparable to benchmark	Better or comparable to benchmark	Benchmark
Cell lysis efficiency for viral vector processing	Better or comparable to benchmark	Better or comparable to benchmark	Benchmark
Endotoxin removal for plasmid manufacturing (with Natrix® Q chromatography membrane) detergent as Wash Buffer Supplement	99.997% removal	99.969% removal	99.997% removal
Emprove® & M-Clarity™	Emprove® Evolve, MQ400, (ISO9001 or equivalent)	Emprove® Expert, MQ500, (IPEC-PQG GMP)	Emprove® Expert, MQ500, (IPEC-PQG GMP)
Composition	30% aqueous solution	Pure substance (100%)	Pure substance (100%)
Available formats	1 L, 2.5 L, 25 L	1 L, 2.5 L, 25 L	1 L, 2.5 L, 25 L, 190 L
Customized Certificate of Analysis options	Available	Available	Available
Toxicology Report	Available	Available	NA

Reference: <sup>1</sup> <https://echa.europa.eu/>

The Life Science business of Merck operates as MilliporeSigma in the U.S. and Canada.

## Choosing Fit-For-Purpose Raw Materials

As a supplier of thousands of raw materials for biologics manufacturing, we know the importance of quality and consistency. Our M-Clarity™ and Emprove® programs are designed to provide customers with the transparency and information needed to select fit-for-purpose raw materials and access to detailed documentation to meet regulatory requirements.

To learn more about our M-Clarity™ attributes, click here or download our brochure

For our Emprove® Program click here or download our brochure



To place an order or receive technical assistance in Europe, please call Customer Service:

France: 0825 045 645

Germany: 069 86798021

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: [MerckMillipore.com/offices](https://www.MerckMillipore.com/offices)

For Technical Service visit: [MerckMillipore.com/techservice](https://www.MerckMillipore.com/techservice)

Merck KGaA

Frankfurter Strasse 250

64293 Darmstadt, Germany

[MerckMillipore.com](https://www.MerckMillipore.com)

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

**Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®**

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q, SAFC, Sigma-Aldrich, Supelco, Mark<sup>1</sup>, Mark<sup>2</sup>, and Mark<sup>3</sup> are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK\_FL12567EN Ver. 1.0 49745 07/2023

