

LEUKOMAX® Leukopaks

Research and GMP Grade

BioIVT offers research grade LEUKOMAX® and GMP grade VivoSTART™ leukopak units from normal and disease state human donors collected at our donor centers under IRB-approved regimens.



With our global network of donor centers, your research doesn't have to wait!

- Same-day, fresh leukopak deliveries within 200 miles of our 11 donor centers in the U.S and U.K.
- Customize your collection with peripheral blood mononuclear cell (PBMC) or immune cell subset isolations
- Our VivoSTART[™] GMP-grade leukopaks are collected in BioIVT-owned, GMP-compliant facilities adhering to 21 CFR 1271. For more details, see our VivoSTART[™] Leukopak Product Guide
- Donor demographic information available, including:
 - Age
 - HLA-type
 - Ethnicity
- Gender
- Smoking status
- Medications & surgical history

Over 600 HLA-typed recallable donors who have donated in the past 6 months, with over 250 being HLA 0201.

Available in Fresh and Cryopreserved Formats

Leukopak Size	Total Number of PBMCs*	Fresh	Cryopreserved	
Quarter Pak	1-2 Billion	HUMANLMX25-0001127	HUMANLMX25-0001124	
Half Pak	3-4 Billion	HUMANLMX50-0001128	HUMANLMX50-0001125	
Full Pak	6-8 Billion	HUMANLMX100-0001129	HUMANLMX100-0001981	
Screening Pak A	500 Million	HUMANLMX500M-0103521	N/A	
Screening Pak B	1 Billion	HUMANLMX1B-0002706	N/A	

Cell Type	Percentage per Leukopak*			
T Cells	57%, range 45-60%			
Monocytes	20%, range 10-30%			
B Cells	10%, range 5-15%			
NK Cells	7.5%, range 5-10%			
Granulocytes	5.2%, range 0.6%-10%			
CD34+ Stem Cells	0.3%, range 0.1-0.5%			

* Values estimated

Your Complete Cell & Gene Therapy Starting Material Provider



Maximize your cell count with mobilized leukopaks

Mobilized LEUKOMAX® leukopak units provide a rich source of CD34+ hematopoietic stem cells compared to traditional leukopaks. Donors are prepped with:

- Granulocyte-colony stimilulating factor (G-CSF), also known as Neuopgen[®], FDA approved
- Plerixafor, also known as Mozobil[®]

Description	Average* WBC Yield	Average * CD34+ Yield	Days					
			1	2	3	4	5	6
Mozobil X1 Collect Day 2	2.5 x 10 ¹⁰	1.2 - 2.5 x 10 ⁸	Mozobil	Collect				
Mozobil X2 Collect Days 2 & 3	5 x 10 ¹⁰	2.5 - 5 x 10 ⁸	Mozobil	Mozobil Collect	Collect			
Neupogen X5 Mozobil X2 Collect Days 5 & 6	9 x 10 ¹⁰	4.5 - 9 x 10 ⁸	Neupogen	Neupogen	Neupogen	Neupogen Mozobil	Neupogen Mozobil	Collect
Neupogen X5 Mozobil X1 Collect Days 5 & 6	8 x 10 ¹⁰	4 - 8 x 10 ⁸	Neupogen	Neupogen	Neupogen	Neupogen Mozobil	Neupogen Collect	Collect
Neupogen X5 Days 5 & 6	7 x 10 ¹⁰	3.5 - 7 x 10 ⁸	Neupogen	Neupogen	Neupogen	Neupogen	Neupogen Collect	Collect

* Values estimated


Disease State Leukopaks

Want to get closer to the real-world starting material? BioIVT offers disease state leukopaks from select indications.

- Physician confirmed diagnosis
- Disease progressions, treatment history and co-morbidities
- Recallable donors
- Quarter and Half Leukopak volumes available

For more information on BioIVT's donor screening, viral testing, and collection processes, visit our leukapheresis collections page





Disease State Leukopak Availability



North America & Asia Pacific T +1 516 483 1196 F +1 516 483 4683 E customerservice@bioivt.com Europe, Middle East & Africa T +44 1444 707333 F +44 (0) 1444 250066 E cseurope@bioivt.com