


FARRAR® REFRIGERATION SOLUTIONS FOR LIFE SCIENCE APPLICATIONS

MODULAR STORAGE



ULC-190

- **Application:** (+2°C to +8°C and -20°C to -80°C operation) redundant precision forced-air convection freeze/store chamber w/ material handling
- **Storage Volume:** 5,380 L



ULC-259

- **Application:** (+2°C to +8°C and -20°C to -80°C operation) redundant precision forced-air convection freeze/store chamber w/ material handling
- **Storage Volume:** 7,334 L



ULC-328

- **Application:** (+2°C to +8°C and -20°C to -80°C operation) redundant precision forced-air convection freeze/store chamber w/ material handling
- **Storage Volume:** 9,288 L

BIOPROCESSING



4000-LC SERIES CONTROLLED RATE CHAMBER

- **Application:** Precision HMI-controlled forced-air convection processing controlled-rate freeze/thaw chamber
- **Storage Volume:** 659.8 L



4000 CONTROLLED RATE CHAMBER

- **Application:** Precision Analog-controlled forced-air convection processing controlled-rate freeze/thaw chamber
- **Storage Volume:** 659.8 L

LABORATORY



FLR105

- **Application:** 4°C Undercounter precision temperature control refrigerator
- **Storage Volume:** 150 L



FLF105

- **Application:** -30°C Undercounter precision temperature control freezer
- **Storage Volume:** 150 L



FLR120

- **Application:** 4°C Upright precision temperature control refrigerator
- **Storage Volume:** 572 L



FLR125

- **Application:** 4°C upright precision temperature control refrigerator
- **Storage Volume:** 714 L



FLR245

- **Application:** 4°C upright precision temperature control refrigerator
- **Storage Volume:** 1271 L



FLF120

- **Application:** -30°C upright precision temperature control freezer
- **Storage Volume:** 572 L



FLF125

- **Application:** -30°C upright precision temperature control freezer
- **Storage Volume:** 714 L



FUF126

- **Application:** -80°C upright precision temperature control freezer
- **Storage Volume:** 736 L

Connect with a cold chain expert today! Discuss your needs, schedule a product demonstration, receive a quote, or request more information.

Gulsah.Metinoz@tranetechnologies.com

CORPORATE HQ

800-B Beaty Street, Davidson, NC 28036

MANUFACTURING HQ

14400 Bergen Boulevard, Noblesville, Indiana 46060

EUROPE HQ

Lenneke Marelaan 6, 1932 Sint-Stevens-Woluwe Belgium