

Condensed Matter

School of Physics and Astronomy

HENRY ····
ROYCE ····
INSTITUTE



UNIVERSITY OF LEEDS

Growth Techniques

IOP Magnetism Winter School Lecture

Dr Philippa Shepley

Condensed Matter


School of Physics and Astronomy

HENRY
ROYCE
INSTITUTE




UNIVERSITY OF LEEDS


Royce Deposition System



Thin film materials




<1 nm to 100 nm




- We combine layers of different materials in “multilayers”.
- Multilayers allow us to study the interfaces between materials.

Royce Deposition System




MBE
(topological
insulators)




PLD
(Complex
oxides)

- Chambers linked under UHV
- Multilayer thin film deposition
- Physics at interfaces



Sputtering
(Metals and
oxides)



MBE
(Organics
and metals)

Outline



UNIVERSITY OF LEEDS


- How do we achieve, maintain and measure vacuum?
 - Vacuum-compatible materials
 - Seals
 - Cleaning
 - Baking
 - Load locks
 - Pumps
 - Gauges
- Thin film deposition processes and techniques
 - Sputtering
 - Pulsed Laser Deposition
 - Molecular Beam Epitaxy
- In situ measurements

Atmosphere to UHV



UNIVERSITY OF LEEDS

- **Vacuum-compatible materials**
 - Use materials with low vapour pressures to limit evaporation into the chamber
 - Different materials are suitable for different levels of vacuum and temperatures
- Stainless steel (chamber walls)
- Oxygen-free copper (gaskets)
- High melting point metals
 - e.g. molybdenum, tantalum, tungsten
- High melting point ceramics
 - e.g. alumina, boron nitride



Vapor pressure vs. temperature



HV
UHV

Atmosphere to UHV



UNIVERSITY OF LEEDS

- **Seals**
 - Prevent gas leaking into the chamber from atmosphere
 - Different seals for different levels of vacuum
- Copper gaskets for UHV seals
- Knife edge on flange cuts into copper as bolts are tightened




Atmosphere to UHV



UNIVERSITY OF LEEDS

- **Cleaning**

- Clean all components and substrates by sonicating in acetone to remove organic residues, then isopropanol to remove the acetone
- Use gloves and tweezers to handle components and substrates



- An ultrasound bath introduces high frequency vibration to help dislodge dirt, grease etc

Pumps





UNIVERSITY OF LEEDS




- Pumps**

- Remove gas molecules from the system
- Different pumps used for different levels of vacuum
- Different pumps for different applications and gases

0.75 Torr
= 1 mbar
= 100 Pa




Turbomolecular pump 10^{-2} - 10^{-10}



Diaphragm pump
 $>10^3$ - 10^{-1}

Cryogenic pump
 10^{-3} - 10^{-9}





Atmosphere to UHV



UNIVERSITY OF LEEDS

- **Baking**

- Heat the entire chamber to $>100^{\circ}\text{C}$ to speed up evaporation of volatile molecules (e.g. H_2O) from the chamber walls
- After baking, individual components can be slowly heated to **outgas**



- **Load locks**

- Using separate air-lock chambers to load and unload samples means we don't regularly vent the whole chamber



Pressure gauges




UNIVERSITY OF LEEDS



- **Gauges**

- Measuring pressure to understand processes
- Measuring pressure to take action to get the right pressure
- Different gauges used for different levels of vacuum


**0.75 Torr
= 1 mbar
= 100 Pa**




Baratron
10⁻¹-10⁻⁴



Ion gauge 10⁻⁶-10⁻¹¹



Pirani gauge 10³-10⁻³





Cold cathode gauge 10⁻²-10⁻⁹

Residual Gas Analysis



UNIVERSITY OF LEEDS

- Mass spectroscopy of residual gases
- Gives partial pressures of different gases by mass
- Can check for leaks by applying He outside chamber monitoring He partial pressure



Pressure range
 10^{-4} - 10^{-11} mbar

Deposition techniques



UNIVERSITY OF LEEDS



- Different techniques are useful for different materials and applications.
- Sputtering
 - Magnetron
 - dc/rf
 - Off axis sputtering
- Pulse laser deposition (PLD)
- MBE
 - Thermal evaporation
 - Electron-beam

Magnetron sputtering



UNIVERSITY OF LEEDS

Sputter chamber





- Metals, oxides, nitrides
- Ideal for polycrystalline or textured films, but epitaxy is possible
- Atoms impinge on substrate with high energy
- Typical growth rates $\text{\AA}/\text{s}$
- Very scalable and common across many industries

Magnetron sputtering



UNIVERSITY OF LEEDS

- A dc voltage between the target and shield ignites an Ar plasma e.g. voltage 340 V (current 0.24 A, power 82 W)
- Permanent magnets behind the target create a magnetic field that confines the plasma into a torus
- Ar ions from the plasma hit the target, ejecting (“sputtering”) atoms and small clusters of atoms
- Sputtered particles are deposited onto the substrate to form a thin film e.g. rate of deposition 1 Å/s
- The rate of deposition and (crystal) structure of the film depend on the gun power, process gas and T-S distance (pressure x distance), the target material and the substrate temperature




Magnetron sputtering




UNIVERSITY OF LEEDS

- Vacuum chamber at 10^{-8} to 10^{-10} mBar
- Load a substrate plate into the load-lock chamber and pump
- Transfer the substrate to the main chamber
- Set the substrate temperature e.g. 25°C
- Fill the main chamber with process gas e.g. Ar ~ 3 μBar



- Apply voltage between the target and shield of the magnetron
- Open the magnetron shutter for a set deposition time
- Deposit a multilayer thin film by using different magnetrons
- Ensure the sample is at room temperature
- Move the sample to the load-lock and out of the system




Sputtering



UNIVERSITY OF LEEDS

Radio frequency (rf) sputtering


- Use an rf voltage to create a plasma
- Avoid build up of charge on insulating targets
- Typically used for insulating targets e.g. oxides



AC electric field
between cathode
and anode
kHz - MHz

Off-axis sputtering

- Capture atoms that are scattered out to the side
- Lower energy atoms impinging on the substrate
- Low energy \rightarrow low damage \rightarrow better epitaxy
- Slower growth rates $\sim 0.01 \text{ \AA/s}$




Pulsed Laser Deposition



UNIVERSITY OF LEEDS



Pulsed laser deposition (PLD)

complex
oxide
thin film
substrate



SrTiO_3
 SrRuO_3
 BiFeO_3

- multi-target stage for multilayer capability



- Good transfer of stoichiometry
- Can be used to grow a range of materials, but commonly used for complex oxides
- Epitaxial films over a small area
- Typical growth rates 0.1 \AA/s

Pulsed Laser Deposition



UNIVERSITY OF LEEDS

- A pulse of laser light transfers energy to the target surface
- The temperature of the target surface increases
- Material from the target surface evaporates or sublimates
- The vaporised (“ablated”) material forms a plasma
- The plasma expands in a “plume”


- There are a lot of processes happening at the target surface and in the plume!

- The properties of the plume depend on the **target material**, **laser fluence** and **process gas pressure**

- Ions and small particles of target material from the plume hit the substrate surface to form a thin film
e.g. rate of deposition 0.1 Å/pulse

- The rate of deposition and (crystal) structure of the film depend on the properties of the plume, T-S distance and the temperature of the substrate

The ablation process



Pulsed Laser Deposition



UNIVERSITY OF LEEDS

- Vacuum chamber at 10^{-8} to 10^{-10} mBar
- Load a substrate plate into the load-lock chamber and pump
- Transfer the substrate to the main chamber
- Fill the main chamber with process gas e.g. O_2 ~ 0.1 mBar
- Heat the substrate to the growth temperature e.g. $700^\circ C$
- Measure the fluence of the laser with an energy meter

$$\text{e.g. fluence} = \frac{\text{pulse energy}}{\text{area on target}} = \frac{30 \text{ mJ}}{0.02 \text{ cm}^2} = 1.5 \text{ J/cm}^2$$


- Fire the laser to pre-ablate the target (clean the target surface)
- The target is scanned during ablation so that we use the whole surface



248 nm Kr-F
excimer laser
pulse

plume

ablated
area






Pulsed Laser Deposition



UNIVERSITY OF LEEDS

- Open the shutter and ablate the target to deposit material onto the substrate
 - e.g. 21000 pulses at a pulse rate of 5Hz
- The substrate sits at the edge of the plume
 - e.g. target-sample distance (T-S) 55 mm
- We can measure the crystal structure of the surface during growth using RHEED
- Cool the sample in a process gas
 - e.g. O₂ at 200mBar
- Transfer the sample to the load-lock and out of the system



Molecular Beam Epitaxy




UNIVERSITY OF LEEDS


Molecular Beam Epitaxy (MBE)

Epitaxial thin film

substrate





Bi_2Se_3



Cu/C_{60} multilayers

- Beam of evaporated/sublimed material




- Wide range of source materials: metals, organics
- Very specialized chambers dedicated to a few materials
- Alloys are made by combining different sources,
- Stoichiometry can be challenging
- Very slow growth rates often $<0.1 \text{ \AA/s}$
- Potential for very good epitaxy and high degree of crystallinity

MBE: effusion cell




UNIVERSITY OF LEEDS

- Thermal evaporator with a crucible made of e.g. boron nitride, tungsten, graphite, quartz
- The material in the crucible is heated using resistive heating elements until a vapour forms
- “effusion” is the process of atoms/molecules as the vapour exiting an aperture
- Cells designed for higher temperatures are needed for materials with low vapour pressures
- Can include multiple heaters e.g. heat the crucible lip to a higher temperature to prevent vapour from condensing there
- Option to introduce reactive gas or plasma for e.g. oxide growth



Cracker cell

- Some materials tend to be evaporated as molecules from standard effusion cells
- Use a heated valve to control vapour output
- Use a heated cracking tube to break up molecules





Schematic sketch of the VCCS building blocks. The PBN-funnel (2) connects the PBN reservoir (1) with the full-PBN valve (3). The evaporant is guided to the cracking stage (5) by a PBN injector tube (4). For better homogeneity there is a beam shaping funnel (6) on the orifice of the VCCS.

MBE: e-beam evaporation




UNIVERSITY OF LEEDS

- Electron beam emitted from a filament
- Beam steered onto the target material with magnetic field
- Electron bombardment heats and melts the target
- Works well for high melting point materials



- Hearth can contain multiple pockets for different target materials

- UHV technique 10^{-8} - 10^{-11} mbar
- Prepare chamber – sources have to be pre-heated to so that unwanted gases can be pumped away
- Heat substrate to outgas
- Might flash-anneal to reconstruct the substrate surface
- Shutter covers the substrate
- Heat sources (effusion/e-beam) until the right flux of material is being emitted
- Open the shutter to start depositing
- RHEED is often used to monitor growth
- Close the substrate shutter
- Cool sources and substrate




In situ measurement



UNIVERSITY OF LEEDS



Flux monitors

- Measuring the partial pressure of material vapour
- Ion gauge filament
- Used for MBE in UHV conditions
- Position at the source focal point
- Retract flux monitor and move substrate down for growth



Quartz crystal microbalance

- Measuring thickness quantity of deposited material
- Apply ac voltage to oscillate quartz
- Measure frequency shift due to deposited material








In situ measurement - RHEED




UNIVERSITY OF LEEDS

Reflection High Energy Electron Diffraction (RHEED)


- Electron diffraction at grazing incidence ($2\theta \sim 1^\circ\text{-}5^\circ$)
 - Highly surface sensitive
- Typical accelerating voltages 10-30 kV
- Electron wavelength $\ll 1 \text{ \AA}$
 - Sensitive to atomic length scales
- Real-time monitoring of layer formation



Si(7x7)
surface
pattern



sub 1ML
pattern gone



Epi-film,
streaks


Summary




- How we achieve, maintain and measure vacuum
 - Cleaning
 - Vacuum-compatible materials
 - Seals
 - Baking
 - Load locks
 - Pumps
 - Gauges
- Thin film deposition processes and techniques
 - Sputtering
 - Pulsed Laser Deposition
 - Molecular Beam Epitaxy
- In situ measurements

Royce Deposition System advice and access:
Dr. Philippa Shepley p.m.shepley@leeds.ac.uk
condensed-matter.leeds.ac.uk

Royce Deposition System




MBE
(topological
insulators)




PLD
(Complex
oxides)

- Chambers linked under UHV
- Multilayer thin film deposition
- Physics at interfaces



Sputtering
(Metals and
oxides)



MBE
(Organics
and metals)

Capabilities



- Thin film deposition with precise control of layer thickness (<1 nm to ~100 nm)
- Substrates sizes of a few mm up to 2 inch diameter wafers
- Wide range of materials (metals, oxides, organics) with best suited deposition techniques (sputtering, pulsed laser deposition, molecular beam epitaxy)
- Temperature control of substrates -100 to 1000°C
- Ultra high vacuum transfer between chambers for clean interfaces between different materials
- RHEED for in-situ characterisation of surfaces

Access for samples via the Royce Access Schemes: www.royce.ac.uk/access-schemes/

Contact me to discuss access: p.m.shepley@leeds.ac.uk

Royce Deposition System

HENRY ····
ROYCE ····
INSTITUTE



UNIVERSITY OF LEEDS

17th-20th March, Imperial College London

ROYCE TRAINING: THIN-FILM DEPOSITION

Applications opening soon! www.royce.ac.uk/events/