

## **ATEC**®

## Breast Biopsy & Excision System



# The simple answer to a flexible biopsy solution.

Our ATEC® Sapphire 200 vacuum-assisted breast biopsy console is a simple, all-in-one platform for all three modalities: ultrasound, MRI, stereotactic 2D and 3D™-image guided biopsy. Designed for radiologists and surgeons who need a simple solution with maximum utility, that considers compassionate care and results in essential elements in breast biopsy.

Designed for compassionate patient care, the system delivers consistent, quality cores in one needle insertion.





## ATEC Breast Biopsy Needle for Ultrasound

A proven solution that enables a fast and compassionate breast biopsy.

Powerful tethered vacuum delivers multiple samples with a single insertion and local anesthetic without interruption.<sup>1</sup>



### ATEC Breast Biopsy Needle for MRI

The pioneer in breast MRI biopsy procedures, designed for compassionate patient care, speed and simplicity.

Delivering the fastest procedure time in MRI by 10 minutes.<sup>2</sup>



### **Eviva® Breast Biopsy Needle** for Stereotactic

Seamless integration with the ATEC Breast Biopsy device, provides optimised access to challenging lesions with biopsy needle options for most breast sizes.

## Technical **specifications**

Electrical Requirements	
Input Line Voltage	220-230VAC
Maximum Power	1150W
Frequency	50/60Hz
Maximum Current	5Amps
Fuse	6Amps, Breaker

ATEC® Sapphire 200 Console		
Construction Dimensions	Caster wheels with 360 degrees rotation 25 x 21 x 37in (63 x 53 x 93cm) WxDxH	
Weight	110lbs (50kg)	
Safe Working Load	22lbs (10kg)	
Power Cord Length	19.8 or 15ft (3 or 5m)	
Vacuum Generated	~28"Hg (71cmHg) at sea level	
Footprint	3400cm	
Cycle Time	4.5 seconds for acquisition	





ATEC Handpiece Assembly



ATEC & Eviva Biopsy Devices

Operating Environmental Conditions		
Ambient Temperature	+59 to + 104°F (+15 to +40°C)	
Relative Humidity Range	30% to 75%	
Atmostpheric Pressure Range	812 hPa to 1014 hPa	

ATEC Footswitch - MRI Conditional	
Size	4 x 6 x 2in (9 x 16 x 5cm) WxLxH
Weight	1lb (0.5kg)
IPX Rating	IPX8
Cable Length	20ft (6m)
Accessories	
ATEC Sapphire 200	Canister with lid
ATEC Biospy Device	Tissue Filter Remote Tissue Filter (RTFA) Stereotactic Adapter Needle Guides Introducer Localization System (ILS)

Tissue Filter

Stereotactic Adapter Needle Guides

ATEC Biopsy Device	
Handpiece Size	1.66 diameter x 10.56in long (4.22 x 26.83cm)
Handpiece Weight	7.26oz (204g)
Needle Lengths	3.52in (9cm) 4.72in (12cm) 5.50in (14cm)
Outer Cannula Diameter	12g: 0.111in (2.82mm) 9g: 0.148in (3.76mm)
Inner Cannula Diameter	12g: 0.084in (2.13mm) 9g: 0.118in (3mm)
Aperture Length	0.787in (20mm) 0.472in (12mm)
Tubing Set Length	US/STX handpiece: 12ft (3.66m)
	MRI handpiece: 20ft (6.10m)
Sterilisation	By radiation
Eviva Biopsy Device	

Eviva Biopsy Device	
Handpiece Size	1.48 x 8.17 x 1.81in (3.76 x 20.75 x 4.6cm) WxLxH
Handpiece Weight	8.7oz (206.6g)
Needle Lengths	3.93in (910cm) 5.11in (13cm)
Outer Cannula Diameter	12g: 0.111in (2.82mm) 9g: 0.148in (3.76mm)
Inner Cannula Diameter	12g: 0.080in (2.03mm) 9g: 0.118in (3mm)
Aperture Length	0.787in (20mm) 0.472in (12mm)
Tubing Set Length	STX handpiece: 12ft (3.66m) MRI handpiece: 20ft (6.10m)
Sterilisation	By radiation





Eviva® Biospy Device

Hologic BV, Da Vincilaan 5, 1930 Zaventem, Belgium.

Breast & Skeletal Health | Hologic.com | euinfo@hologic.com



References: 1. Schrading S, Simon B, Braun M et al. MRI-Guided Breast Biopsy: Influence of Choice of Vacuum Biopsy System on the Mode of Biopsy of MRI-Only Suspicious Breast Lesions. AJR. 2010. 194:1650-1657. 2. Compared to Mammotome and Vacor. M. Scarth, W. Teh; Harrow/UK. MR-guided vacuum-assisted core-needle breast biopsy: Comparison of three vacuum-assisted biopsy devices. ECR Congress 2009. Scientific Paper, e-Poster: B-295

DS-10219-EUR-EN REV. 001 ©2023 Hologic Inc. All rights reserved. Hologic, The Science of Sure, ATEC, Eviva and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to euinfo@hologic.com.