



E.An-II Plus

3" High Resolution Full
Angle 16:9 Screen
Rechargeable Battery
Water Proof (IP 44)
Nishikawa Blade
Camera & Video
Anti Fog







FIGURE 4–8 Nishikawa Laryngoscope blade. This blade has an S-shaped spatula and bifid tip, designed to manage the bulk of a large tongue. (From Nishikawa K, Yamada K, Sakamoto A. A new curved laryngoscope blade for routine and difficult tracheal intubation. *Anesth Analg*. 2008;107:1248–1252 with permission)



Quote from ATLAS OF
AIRWAY MANAGEMENT
TECHNIQUES AND TOOLS

P 33

Features of Nishikawa Blade

- S-shaped Spatula : Designed to manage the enlarged tongue base.
- Bifid Tip: Fits better into vallecula space & lift epiglottis more effectively.



Ordering Information

Product Shot	Product Model	Size	Package	Max. Tracheal Tube Size / Best Patient Weight
	E.AN-IIL E.AN-IIS	Standard Adult Standard Pediatric	1	40-120kg 1.5-50kg
	E.AN-IILP E.AN-IISP	Premium Adult Premium Pediatric	1	40-120kg 1.5-50kg
REGULAR BLADES				
	ZH-L	114mm/16mm	200/case, 50/box	50-120kg
	ZH-M	106mm/15mm	200/case, 50/box	40-100kg
	ZH-S	83mm/14mm	200/case, 50/box	4-50kg
	ZH-N	79mm/9mm	200/case, 50/box	1.5-5kg
NISHIKAWA BLADES				
	ZH-LS	114mm/16mm	200/case, 50/box	50-120kg
	ZH-MS	106mm/15mm	200/case, 50/box	40-100kg
	ZH-SS	83mm/11mm	200/case, 50/box	4-50kg
CHANNEL BLADES				
	ZH-LG	114mm/20mm	100/case, 25/box	Φ 8.5 50-120kg
	ZH-MG	106mm/19mm	100/case, 25/box	Φ 7.5 40-100kg

Medan Video Laryngoscope + Nishikawa Blade

The “BEST” Option for Routine
and Difficult Tracheal Intubation



TEL : +86 022-86839935 +86 4001-1212-10 FAX : +86 022-83717259
E-MAIL : tianjinmedan@yahoo.com tianjinmedan@163.com
ADD : Building 35-1 U-Valley Xinke Garden , Beichen
District, Tianjin, China P.C. 300401 <http://www.tianjinmedan.com>

Tianjin Medan Crop.

medan

Local Distributor: