


# Faxitron® Path Specimen Radiography System

## Work Efficiently With More Imaging Area.

**Introducing the Faxitron® Path Specimen Radiography System from Hologic.**

With a large field of view and high resolution detector, the Faxitron® Path system is designed to deliver greater efficiency and exceptional patient care in the Pathology lab.



### **Faxitron® Path System Features:**

#### **High-Resolution Imaging**

The Faxitron® Path system has a magnification shelf with field of view guides and automatic position detection that provides up to 6X geometric magnification, which can locate the smallest calcifications in tissue samples to ensure accurate diagnoses.

#### **Wide Energy Range**

The system can image both thin slices and large specimens, and proprietary algorithms enable optimal image exposure and clarity with Automatic Exposure Control (AEC).

#### **Improved Efficiency**

Routine access to high-resolution imaging with the Faxitron® Path system results in faster final report generation. The software also has the ability to send multiple annotated images to PACS at the touch of a button.

#### **Optimal Imaging Quality**

The Faxitron® Path system allows the user to optimally image anything from breast tissue slices and intact mastectomies to bone and fetal remains. The system's large detector has a high resolution for its size and accommodates a broad range of samples.

#### **Easy to Use**

No additional training or specialised x-ray requirements are needed to operate the Faxitron® Path system, and the system plugs into any standard A/C outlet.

#### **Convenient Interface**

The continuous surface keyboard allows for ease of use and the system has the ability to save images in a variety of file formats.

#### **Optimized Design**

The Faxitron® Path system allows for optimal placement at the point-of-care. Its stainless steel chamber was designed with no right angles to allow for easy cleaning.

**Faxitron® Path**  
Specimen Radiography System

# A Comprehensive Family of Solutions

The Faxitron® Path system fits right into Hologic's full range of specimen radiography solutions. Hologic has a complete set of specimen imaging products for Radiology, Surgery and Pathology to help you improve efficiency, optimise the patient experience and deliver the best outcomes possible.

Talk to your representative about the solutions that best suit your needs, ensuring your patients have the best possible care when they need it most.

To learn more about Hologic's solutions for minimally invasive surgery, please visit: [hologicbreastsurgery.com](http://hologicbreastsurgery.com)

**Faxitron® Core  
Specimen Radiography System**



**Faxitron® Trident® HD  
Specimen Radiography System**



**Faxitron® Path+  
Specimen Radiography System**



Faxitron® Core system also known as CoreVision® system.  
Faxitron® Path+ system also known as PathVision® XL system.

[www.hologic.com](http://www.hologic.com) | [euinfo@hologic.com](mailto:euinfo@hologic.com) | +32.2.711.4680

**Hologic BV**  
Da Vincilaan 5  
Building Caprese  
1930 Zaventem  
Belgium  
Fax: +32.2.725.2087

MISC-05519-EUR-EN Rev002 (08/20) Hologic Inc. ©2020 All rights reserved. Hologic, The Science of Sure, CoreVision, Faxitron, PathVision, Trident HD, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshow, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.